



# COLORident



Periodontal curettes and scalers with ceramic color ring for better identification, biocompatible (handle 10m)  
Parodontalküretten und Scaler mit Farbmarkierungsring aus Keramik für eine bessere Identifikation, biokompatibel (Griff 10mm)

-  No. 1/2 and 5/6 front tooth region facial and oral  
Nr. 1/2 und 5/6 Frontzahnbereich fazial und oral
-  No. 7/8 side tooth region facial and oral  
Nr. 7/8 Seitenzahnbereich fazial und oral
-  No. 11/12 or 15/16 mesial premolars and molars  
Nr. 11/12 oder 15/16 mesial Prämolaren und Molaren
-  No. 13/14 distal premolars and molars  
and also 17/18 molars  
Nr. 13/14 distal Prämolaren und Molaren  
und ergänzend 17/18 Molaren



COLORident GRACEY 1-2  
lingual / buccal  
404/1-2-ED ©



COLORident GRACEY 5-6  
lingual / buccal  
404/5-6-ED ©



COLORident GRACEY 11-12  
mesial  
404/11-12-ED ©




COLORident GRACEY 13-14  
distal  
404/13-14-ED ©



COLORident GRACEY 15-16  
mesial  
404/15-16-ED ©

## DEEP SCALE GRACEYS

NEW: LEIBINGER Deep Scale Graceys, especially for the use in deep, narrow pockets and for root planing. On these special curettes, the 1st shaft is 3mm longer and the scaling tip is 50% shorter than regular Graceys.



## DEEP SCALE GRACEYS

NEU: LEIBINGER Deep Scale Graceys sind besonders für den Einsatz in tiefen, engen Taschen und für die Wurzelglättung geeignet. Im Vergleich zu den Standard Graceys ist ihr 1er Schaft 3mm länger und gleichzeitig ist das Arbeitende um 50% verkürzt.



COLORident GRACEY 1-2M DEEP SCALE  
lingual / buccal  
404/1-2M-ED ©



COLORident GRACEY 5-6M DEEP SCALE  
lingual / buccal  
404/5-6M-ED ©




COLORident GRACEY 7-8M DEEP SCALE  
lingual / buccal  
404/7-8M-ED ©



COLORident GRACEY 11-12M DEEP SCALE  
lingual / buccal  
404/11-12M-ED ©



COLORident GRACEY 13-14M DEEP SCALE  
lingual / buccal  
404/13-14M-ED ©



COLORident GRACEY 15-16M DEEP SCALE  
lingual / buccal  
404/15-16M-ED ©



COLORident M23  
universally applicable / für alle Bereiche einsetzbar  
404/M23-ED ©



COLORident Scaler 204SD  
Sickle Scaler / Sichel Scaler  
404/204SD-ED ©



COLORident 204S  
Sickle Scaler / Sichel Scaler  
404/204S-ED ©



COLORident HYGENIST  
Sickle Scaler / Sichel Scaler  
404/H6H7-ED ©



COLORident COLUMBIA 2RL  
anterior application / für den Frontzahnbereich  
404/2RL-ED ©



COLORident COLUMBIA 4RL  
posterior application / für den posterioren Bereich  
404/4RL-ED ©



COLORident LANGER 1-2  
mandibular posterior / Unterkiefer posterior  
404/L1-2-ED ©



COLORident LANGER 3-4  
maxillary posterior / Oberkiefer posterior  
404/L3-4-ED ©



COLORident L56  
maxillary and mandibular anterior / Oberkiefer und Unterkiefer anterior  
404/L5-6-ED ©

SHARPCut



Universal-  
Scaler H6-H7

401/H6-H7-ED-TIN  
©



Universal-  
Scaler 204-S

401/204-S-ED-TIN  
©



Universal-  
Scaler 204SD

401/204SD-ED-TIN  
©



Universal-  
Scaler 204-2RL

401/204-2RL-ED-TIN  
©



Universal-  
Scaler 204-4RL

401/204-4RL-ED-TIN  
©



GRACEY 1-2M  
DEEP SCALE  
lingual / buccal  
400/1-2M-ED-TIN  
©



GRACEY 5-6M  
DEEP SCALE  
lingual / buccal  
400/5-6M-ED-TIN  
©



GRACEY 7-8M  
DEEP SCALE  
lingual / buccal  
400/7-8M-ED-TIN  
©



GRACEY 11-12M  
DEEP SCALE  
mesial  
400/11-12M-ED-TIN  
©

- special hardness through special coating
- long lasting sharpness
- less wear out
- spezielle Hartstoffbeschichtung
- lang anhaltende Schärfe
- weniger Verschleiß

ED = ERGODESIGN Handle  
ED = ERGODESIGN Griff



Scalor M23

401/M23-ED-TIN  
◎



GRACEY 1-2  
lingual / buccal  
400/1-2ED-TIN  
◎



GRACEY 3-4  
lingual / buccal  
400/3-4ED-TIN  
◎



GRACEY 5-6  
lingual / buccal  
400/5-6ED-TIN  
◎



GRACEY 13-14M  
DEEP SCALE  
distal  
400/13-14M-ED-TIN  
◎



GRACEY 7-8  
lingual / buccal  
400/7-8ED-TIN  
◎



GRACEY 11-12  
mesial  
400/11-12ED-TIN  
◎



GRACEY 13-14  
distal  
400/13-14ED-TIN  
◎



GRACEY 15-16  
distal  
400/15-16ED-TIN  
◎